

*Application  
for  
United States Letters Patent*

To all whom it may concern:

Be it known that

**Paul E. Harris and Charles Hesdorffer**

has invented certain new and useful improvements in

**GROWTH OF HUMAN DENDRITIC CELLS  
FOR CANCER IMMUNOTHERAPY  
IN CLOSED SYSTEM USING MICROCARRIER BEADS**

of which the following is a full clear and exact description.



antitumor immunotherapy, 4 Cytokines Cell Mol. Ther. 265-73 (1998)}. The open system is labor intensive and poses an increased risk of microbial contamination to the expanded product, the patient and the technician.

An alternative to the open flask is a closed system for culturing populations of monocyte enriched peripheral blood mononuclear cells using flexible gas permeable cell culture bags and sterile connecting devices {A. Glaser et al., *Collection of mononuclear cells in the Spectra for the generation of dendritic cells.* 39 Transfusion 661-62 (1999)}. Growing human DC in plastic bags, even under clinical grade and using good manufacturing practices, have poor yields because the surface of the bags is suboptimal.

#### SUMMARY OF THE INVENTION

The present invention, in accordance with one embodiment, provides a method of reproducibly generating dendritic cells, comprising the steps of (a) loading blood mononuclear cells into a cell culture container containing microcarrier beads therein, (b) incubating for a predetermined time period tissue culture comprising the cells loaded in the container in step (a), and (c) separating nonadherent cells and cells adhered to the beads.

The present invention, in accordance with another embodiment, provides a method of reproducibly generating dendritic cells, comprising the steps of (a) loading microcarrier beads into a cell culture container, (b) loading blood mononuclear cells into the container, (c) incubating for a predetermined time period tissue culture comprising the mononuclear cells loaded in the container in step (b), and (d) separating nonadherent cells and cells adhered to the beads.







unit of culture surface area in closed gas permeable cell culture bags is less than the yield in open flask systems. To improve the yields of DCs in a closed system, styrene copolymer beads (e.g., 90-500 micron diameter, density  $\geq$  1.04 g/cm<sup>3</sup>, SoloHill Engineering, Inc., Ann Arbor, MI) are introduced into the bags, in accordance with one embodiment of the present invention, to increase the available surface area (e.g., by 380 cm<sup>2</sup>) and supply a surface area similar to that found in the flasks {see, e.g., M. Kiremitci et al., *Cell adhesion to the surfaces of polymeric beads*, 18 Biomater. Artif. Cells Artif. Organs 599-603 (1990)}. Also, the beads have a density that allows them to sink/settle in due course.

In accordance with one (FIGS. 3 through 5) of many possible embodiments, 1 gram of gamma radiation sterilized beads and 10 x 10<sup>8</sup> total cells/bag of MNC product are diluted in 100 mls AIM-V (e.g., from GIBCO, Grand Island, New York) and loaded into gas permeable cell culture bags (e.g., Lifecell X-fold Cell Culture Containers PL2417, 180 cm<sup>2</sup>, Nexell Therapeutics, Irvine, CA), under a biological safety cabinet (step 402).

The tissue culture bags then are incubated (step 403), for example, in a humidified 37°C, 5% CO<sub>2</sub> atmosphere for approximately four hours. After four hours the contents of the bag are gently resuspended (step 404), the beads are allowed to settle for 5 minutes at 1 x g (step 405), and the bag is clamped 1 cm above the settled beads (step 406). The supernatant then is expressed off (step 407) using a transfer bag and a sterile connecting device (e.g., from Terumo Corp., Phoenix, AZ). This procedure (steps 402-407) may be repeated three times with 50 mls AIM-V media. As control for adherence, a sample of the expressed cells may be immunophenotyped for monocyte markers (e.g., CD14 and CD11c). Adherence of MNC to the bag and bead surface may





The quantities of the cells produced are acceptable for adaptive transfer strategies. Current tumor antigen vaccine protocols require approximately  $10^7$  to  $10^8$  total DCs. Using this closed system of culture, a sufficient number of DCs can be harvested for a complete course of therapy using a single 10-liter MNC apheresis and an average of five culture bags.

Another embodiment will now be described with reference to FIGS. 6A through 6D.

A dendritic cell culture medium is prepared (step 701) by combining AIM V media (e.g., BB-MF 2557, Life Technologies, Grand Island, NY), rh-IL-4 (e.g., 1000 U/ml, GLP grade, Sigma Aldrich, St. Louis, MO) and rh-GM-CSF (e.g., 25 ng/ml, Therapeutic grade, Immunex, Seattle, WA).

Polystyrene copolymer beads (e.g., 250 micron diameter, density =  $1.07 \text{ g/cm}^3$ , BB-MF 3094, Solohill Engineering, Inc., Ann Arbor, MI) are obtained and prepared (step 702) for use. For example, the beads may be suspended in a phosphate buffered saline (e.g., EDR9865, therapeutic grade, Nexell, CA) [100 gms beads/200 mls saline], placed in an autoclavable glass bottle and capped, and sterilized in an autoclave (e.g., 20 lbs/sq.in. at  $121^\circ\text{C}$  for 1 hour with slow exhaust). The container is sealed and then transferred to a biological safety cabinet. A 1 ml aliquot is removed, placed in the transport tube, and tested for sterility (e.g., Bioscreen Testing Services, Inc., Torrance, CA). The polystyrene beads are used if no bacterial growth is detected.

Also, a sterile peptide (e.g., HER-2, or another peptide antigen specific to another target tumor) solution is prepared (step 703), using for example synthetic peptide (e.g., GLP grade, American Peptide Company) and phosphate

buffered saline. For example, HER-2 synthetic peptide powder is dissolved in saline at a concentration of 200 ug/ml (20x) and sterile filtered through 0.2 micron nylon membrane. The solution is aliquoted in sterile 10 ml vials and stored.

A cryoprotectant agent (e.g., DMSO, USP grade, Gaylord Chemical Corporation, Slidell, LA) is obtained and tested for sterility (step 704).

Apheresis products are transferred to a transfer bag (step 705). Samples of the apheresis products are run through quality control (e.g., hematology analyzer, Trypan blue viability, CD45/14 immunophenotype) [step 706]. If quality control is passed, apheresis products (e.g.,  $10 \times 10^8$  mononuclear cells/bag x four to five bags) are transferred (step 707) from the transfer bag using a sterile connecting device (e.g., Lifecell transfer set, Nexell Therapeutics, Irvine, CA) to gas permeable tissue culture bags (e.g., therapeutic grade, Lifecell X-fold Cell Culture Containers PL2417, 180 cm<sup>2</sup>, Nexell). The beads that pass quality control (step 702) also are inserted in the bags.

The tissue culture bags then are incubated (step 708), for example, in a humidified 37°C, 5% CO<sub>2</sub> atmosphere for approximately four hours. At the midpoint, the bag is flipped from one side to the other. After four hours, the tissue culture is washed three times with AIM-V to remove nonadherent lymphocytes, platelets, grans, RBC, etc. (step 709). The wash includes transfer of the AIM-V media and expressing off the supernatant while leaving the beads in the bag.

Next, the dendritic cell culture medium (prepared in step 701) is transferred via a sterile process to the tissue culture bags (step 710). The bags are incubated again, for

example, in a humidified dedicated 37°C 5% CO<sub>2</sub> incubator (step 711) for 5 to 7 days. At day 4, samples of cell suspension are removed for quality control (step 712). On day 7, the tissue culture bags are moved from the incubator to a biological safety cabinet (step 713). The bags are cooled to room temperature, and the contents are gently agitated for five minutes (step 714). The bags are suspended in an upright position to allow the beads to settle for 5 minutes at 1 x g (step 715), and the bag is clamped above the settled beads (step 716). The cell culture suspension is harvested (step 717) by expression into transfer bags (e.g., Stericell bags, Nexell Therapeutics) using a sterile connecting device (e.g, from Terumo Corp., Phoenix, AZ).

Samples are removed from the transfer bags and run through quality control (step 718). For example, if no beads are present and viability is greater than 95%, then the samples are passed to immunophenotype by flow cytometry. If quality control is passed, cells (e.g., approximately 50 x 10<sup>6</sup> cells or any range, such as all of the cells) can be transferred (step 719) to a second bag for cryopreservation and immunological function controls (e.g., seven day proliferation assay using harvested cells as stimulators for lymphocytes from three different individuals).

HER-2 peptide solution is added to the transfer bag (step 720) for peptide loading onto the HLA Class I of the DC (final concentration is 10 ug/ml). The transfer bag is incubated (step 721) overnight (e.g., 4 to 12 hours at 37°C, 5% CO<sub>2</sub> in a humidified atmosphere of a dedicated incubator). Samples of the peptide loaded DC are tested for mycoplasma (step 722). If the test results are negative for mycoplasma, the peptide loaded DCs are washed three times with therapeutic grade phosphate buffer saline (step 723). In a preclinical phase, samples of the suspension may be

removed for quality control analysis, such as for endotoxin (e.g., USP LAL), fluoride (ion specific electrode) and residual organic solvent (GC-MS).

5 Injection formulation is prepared by resuspending the washed DCs (step 724) at a concentration less than  $10^7$  cells/ml (e.g., 3, 6, 9 or  $12 \times 10^6$  cells/ml) in saline supplemented with 5% autologous serum obtained the same day. Samples are removed for quality control (step 725),  
10 such as Gram stain. If quality control is passed, the injection formulation is cleared for administration and injected within four hours of preparation (step 726).

The remaining cells are cryopreserved. First, the peptide loaded dendritic cells (also DC without peptide) are  
15 suspended (step 731) in a solution of therapeutic grade saline supplemented with 5% autologous serum ( $5 \times 10^6$ /ml). Cryoprotectant agent is added to a final concentration of 10% and placed in sterile NUNC vials (5ml) [step 732]. The  
20 cells are placed in a methanol bath at  $-70^\circ\text{C}$  overnight (step 733), then placed in vapor phase liquid  $\text{N}_2$  storage until use (step 734).

After two to three days of storage, a vial of peptide  
25 loaded DC is retrieved (step 735) from the liquid  $\text{N}_2$  storage for quality control testing. Vials are thawed at  $37^\circ\text{C}$  in a biological safety cabinet (step 736). The cells are washed with AIM-V to remove cryoprotectant (step 737).

30 Aliquots are removed for the following assays:

Viability	(>70%)
Sterility USP	(No growth)
Mycoplasma by PCR	(negative)
MLC test (7 day)	(Stimulate proliferative response greater than 3 X BACKGROUND at responder to stimulator ratio of

35

10 to 1)

Endotoxin (USP LAL) <0.06 EU/ml

5 If the quality control test (step 738) is passed, the  
cryopreserved DC are released for thawing. Vials of  
cryopreserved peptide loaded DC are thawed (step 739) at  
37°C in the biological safety cabinet. Cells are washed  
with therapeutic grade saline three times to remove  
cryoprotectant agent (step 740). Washed DC are resuspended  
10 (step 741) in saline supplemented with 5% autologous serum  
obtained the same day. Samples are removed for quality  
control testing (step 742), e.g., viability staining and  
Gram stain. If viability is greater than 70%, then the  
cells are passed to adjust cell concentration to less than  
15  $10^7$  viable cells/ml (e.g., 6, 9 or  $12 \times 10^6$  viable cells/ml)  
x 1 ml. Next, a Gram stain is applied. If the Gram stain  
is passed, the injection formulation is cleared for  
administration and injected within four hours of  
preparation (step 743).

20 The use of polystyrene beads in closed containers has a  
number of advantages over the use of open flasks for  
reproducibly generating dendritic cells, including  
sterility, risks of exposure for workers, higher yield,  
25 etc. Because of these factors, a closed system, e.g.,  
flexible gas permeable plastic tissue culture bags, is  
preferred over the open flasks. The bags alone, however,  
do not provide an ideal surface for the attachment of DC  
precursor cells (monocytes). The introduction of selected  
30 polystyrene beads into the bags provides a surface that the  
monocytes easily adhere to. Once the monocytes have  
matured into DC, their adherence to the polystyrene  
surface, provided by the beads, is significantly reduced.  
At the end of the culture period, DCs no longer adhere to  
35 the beads and are harvested in the supernatant.

09726883 "113000

The beads are selected based, in part, on their size. Since more surface area is desirable, smaller beads in a larger quantity is preferred to larger beads in a smaller quantity.

Also, the specific gravity of the beads which allows them to settle after a period also contributes to their utility in the above-described methods. Since the monocytes adhered to the beads settle with the beads and thereby separate from the undesired cells (e.g., lymphocytes, platelets, etc.) which are removed by expressing off the supernatant.

Quality control is applied in the methods described above to comply with good manufacturing practices criteria.

While the invention has been described with reference to specific embodiments, the description is not meant to be construed in a limiting sense. It should be understood that the invention is not limited to the precise embodiments described herein.

For example, while the closed system in the embodiments described above use a cell culture bag, other cell culture vessels may be used for the closed container. The container need not have any particular shape. It is preferable, but not essential, to have more than one port on the container to facilitate the transfer of materials in and out of the container. It is important, however, that the container is gas permeable to, for example, O<sub>2</sub> and CO<sub>2</sub>. Also, it is important to maintain a ratio of (beads and container) surface area to container volume that allows the container to hold enough media to support the culture period, so that culture only needs to be fed once, rather than repeatedly.

